

Clinispin horizon 6, 12 & 24 Tube Centrifuges

Versatile centrifuges with 3 default preset cycles
for chemistry, coagulation and urine.



**RAPID HORIZONTAL
SEPARATION.**

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

Clinispin horizon 6, 12 & 24 Tube Centrifuges

Versatile centrifuges with 3 default preset cycles
for chemistry, coagulation and urine.

FEATURES

- **Easy To Use** - simple two-button operation to select a preset cycle or lock to a single cycle
 - 3 preset cycles
 - 6, 12 and 24 tube capacity
 - Erasable label for custom cycles
- **Lid Lights** - flash at the end of each cycle & continue flashing until the lid is opened
- **Automatic Lid Lock** - will not run with lid open & lid cannot be opened while samples are spinning
- **Cool Operation** - unique design of the convection based cooling system protects the samples even under the heaviest workloads
- **Clear Shatterproof Lid** - for easy monitoring
- **Brushless DC Motor** - no routine maintenance
- **Reinforced Guard Bowl** - maximum safety & durability
- **2 Stage Horizontal Rotor** - patented swing-out design produces ideal horizontally separated samples



Clinispin horizon 6
- ideal where space
is at a premium and
consistency is key



Clinispin horizon 12
- ideal for midsize
routine labs with a
repeatability emphasis



Clinispin horizon 24
- ideal for sites with
a larger workload

SPECIFICATIONS

	horizon 6	horizon 12	horizon 24
Tube Capacity	6 tubes – 3-10 mL	12 tubes – 3-10 mL	24 tubes – 3-10 mL 12 tubes – 3-15 mL 6 tubes – 50 mL
Dimensions (H x W x D)	36 x 30 x 23 cm (14" x 12" x 9")	38 x 33 x 23 cm (15" x 13" x 9")	43 x 38 x 23 cm (17" x 15" x 9")
Weight	5.4 kg (12 lbs)	15 kg (34 lbs)	17 kg (39 lbs)
Sound Level	64 dB A	56 dB A	59 dB A
Environmental Range	16 – 32°C	16 – 32°C	16 – 32°C
Voltage	95 – 253 VAC	95 – 253 VAC	95 – 253 VAC
Frequency	50/60 Hz	50/60 Hz	50/60 Hz
Power Requirement	220 Watts	415 Watts	415 Watts
Centrifuge Motor	½ H.P. Brushless	½ H.P. Brushless	½ H.P. Brushless
Max. G-Force	2000 xg	2000 xg	2000 xg
Max. Speed	3800 RPM	3700 RPM	3400 RPM
Cycle Time	1 – 30 minutes (± 2%)	1 – 30 minutes (± 2%)	1 – 30 minutes (± 2%)



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**