

# Clinispin™

Vet Combi

## Multi Purpose Veterinary Centrifuge



**VERSATILE.**  
**COMBINATION ROTOR.**

**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.



# Multi Purpose Veterinary Centrifuge

The Clinispin Vet Combi is a high speed centrifuge equipped with a maintenance free drive, digital display and simple interface for silent and efficient operation. The combination Rotor provides efficient separation of blood and urine samples and accurate Haematocrit determination.

## FEATURES

- 3 pre-defined presets for easy and rapid operation at one press of a button:
  - Serum / Plasma (12,000 RPM for 3 minutes)
  - HCT / PCV (12,000 RPM for 5 minutes)
  - Urine (1,500 RPM for 2 minutes)

Pre-sets can be varied for each centrifugation if required

- Rotor allows for quick and easy sample loading
- Maintenance free drive
- Digital display
- Easy to use
- Quiet operation
- Displays RPM or RCF (g force)



## SPECIFICATIONS

**Motor Type:** Brushless DC Motor

**Maximum Speed:** 12,000 RPM

**Maximum RCF:** 9500xg

**Rotor Volume:** 10 x 2ml tubes and  
10 x 75mm capillary tubes

**Programs:** Serum/Plasma - 12,000 RPM for 3 minutes

HCT/PCV - 12,000 RPM for 5 minutes

Urine - 1,500 RPM for 2 minutes

**Timer:** 1-99 minutes

**Power:** AC220-230V 50/60Hz

**Unit Fuse:** 3A 230V

**Plug Fuse:** 5A

**Weight:** 9kg

**Dimensions:** 300 (L) x 270 (W) x 190 (H) mm

**Permitted Temperature Range:** 5°C - 40°C

**Relative Humidity:** <80%

**Noise (at Maximum RPM):** <65dB(A)



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT  
LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**