

Critical Care / Blood Gas Analyser



- Portable
- Room temperature card storage
- Veterinary specific
- WiFi, Bluetooth connectivity



**RAPID RESULTS,
ANYTIME, ANYWHERE.**



Advanced
Partner

Laboratory Diagnostics
Point of Care Solutions

SIEMENS
Healthineers

epoc Critical Care / Blood Gas Analyser

FEATURES

- Accurate blood gas, electrolyte and critical care biochemistry results in 3 minutes
- Veterinary species specific software
- Room temperature card storage
- 92µl sample, fresh whole blood (arterial, venous or capillary)
- Completely portable with rechargeable battery & bluetooth printer
- Ready to use anytime, anywhere - in-practice, theatre, patient side or equine yard
- Easy to use with touch screen technology
- Wi-Fi and Bluetooth enabled - easy transfer of data
- EDM Lite software and large data storage
- Automatic self calibration and quality checks
- Low cost analyser and cartridges
- Single use, self calibrated test with bar-coded **SmartCard™** technology for improved safety and total quality assurance
- Ability to input all patient data, species, temperature and ventilator settings - analyser stores 2,000 patient results

TEST CARD CONTAINS ALL PARAMETERS

Measured Analytes:

pH, pCO₂, pO₂, Na⁺, K⁺, Cl⁻, Ca⁺⁺, Glu, Lac, HCT, BUN, Creatinine, tCO₂

Calculated Values:

HCO₃, BE (ecf), BE (b), sO₂, Hgb, AGapK



3 Simple Steps to Accurate & Immediate Results



State of the Art Blood Analysis

Advanced
Partner

Laboratory Diagnostics
Point of Care Solutions

SIEMENS
Healthineers



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

WOODLEY
EQUIPMENT COMPANY LTD.