



# Veterinary Glucose Meter

Calibrated specifically for dogs, cats and horses



**SIMPLE, ACCURATE,  
GLUCOSE MONITORING.  
LIFETIME WARRANTY.**

**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.

## KEY FEATURES

- FAD-GDH enzyme test strip
- Strip eject
- Single code strip for all species
- Calibrated specifically for Dogs, Cats & Horses
- Biosensor auto code technology
- Accurate results
- Easy to use
- Vets can train owners to home test
- Small sample (0.7µl) of whole blood required
- Large LCD screen
- 5 seconds test time
- Readable range 1.1 - 33.3 mmol/L
- 500 result automatic memory
- Adjustable alert range alarms
- HCT range 20 - 60%
- Lifetime warranty

**Calibrated specifically for dogs, cats and horses**

## SPECIES SPECIFIC CALIBRATION

- The g-Pet PLUS Veterinary Glucose Meter is calibrated specifically for use in Dogs, Cats and Horses.
- Each g-Pet PLUS meter contains all relevant calibration information to ensure an accurate result is reported for each species.
- Vets no longer need to persist with inaccurate human meters as there is now a low cost, veterinary specific alternative with a lifetime warranty.

## KIT CONTENTS

- g-Pet PLUS Meter
- g-Pet PLUS Test Strips (10)
- Control Solution
- Lancing Device
- Lancets (10)
- AAA Batteries (2)
- Carry Case
- Quick Reference Guide



**Testing Made Easy!**

Insert the test strip  
(Meter will turn on automatically)



Apply Blood



Results in 5 seconds



+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT  
LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**