



Arterial Blood Gas Syringes

Available in 1ml and 3ml sizes, with or without needle



SIMPLE, EFFICIENT
BLOOD SAMPLE COLLECTION.

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.



Arterial Blood Gas Syringes

Available in 1ml and 3ml sizes, with or without needle

FEATURES

- Sterile, single-use device for the collection of whole blood samples
- Luer cap to seal the sample
- Luer lock connection for secure needle attachment
- Breathable plug - double layer porous stone structure, auto-exhaust
- Anticoagulant - ultra-precision sprayed solid lithium heparin, does not dilute blood samples - balanced by calcium ions to ensure accuracy of test results
- Preset system for simple and efficient automatic sampling
- Sharps protection guard on versions with needle to help prevent needlestick injury



SPECIFICATIONS



Validated for use with
epoc Critical Care/
Blood Gas Analyser



Product Code Description

WD0230	1ml without needle, 30 I.U. Calcium ions balanced lithium heparin
WD0231	3ml without needle, 80 I.U. Calcium ions balanced lithium heparin
WD0232	1ml with needle 23G x 25mm with protection, 30 I.U. Calcium ions balanced lithium heparin
WD0233	3ml with needle 23G x 25mm with protection, 80 I.U. Calcium ions balanced lithium heparin

Supplied in Packs of 50

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**