

InSight Hgb / HCT Meter

The InSight Hgb / HCT Meter measures total haemoglobin and haematocrit in capillary and venous whole blood.

- Handheld meter
- Results in less than 12 seconds
- Easy operation
- Up to 250 results stored

Quick Results

The Hgb / HCT Meter displays the total haemoglobin concentration in grams/per decilitre (g/dL) in 12 seconds.

Storage

Up to 250 results are automatically stored with a date and time recorded.

Easy Operation

Simply place a drop of whole blood onto the test strip and the meter's optical detector will automatically measure the change in membrane reflectance.



Three Step Operation

1. Insert strip
2. Apply Blood Drop
3. Read results in 30 seconds

Visit www.woodleylabdiagnostics.com for full product range

InSight Hgb / HCT Meter

Specification

Parameter: Haemoglobin and Haematocrit in fresh whole blood

Measuring Principle: Optical reflectance

Strip: InSight Hgb/HCT Test Strip

Measuring time: Less than 12 seconds

Haematocrit (HCT): Display in percentage

Hgb Measuring range: 4.0g/dL~24.0g/dL

Display: LCD, test results display in SI

Storage Function: Save the strip code, memorise and update 250 test results automatically

Weight: <58g (Battery included)

Operating Storage: 15~30°C ; RH≤80% non-condensing

Non-operating Storage: -20~55°C ; RH≤90% non-condensing



Kit Contents

- InSight Hgb / HCT Meter
- InSight Hgb / HCT Test Strips (10)
- InSight Hgb / HCT Code Card
- Lancets (10)
- Lancing Device
- Instruction Manual
- CR2032 Batteries (2)
- Carry Case

Visit www.woodleylabdiagnostics.com for full product range