

# InSight<sup>®</sup>

## Urinalysis

MS Reader & Urine Test Strips



**QUICK, EASY,**  
**ACCURATE RESULTS.**



**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.

# InSight®

## Urinalysis

Urinalysis is an important tool in disease detection, as well as monitoring and screening animal health. Abnormalities can be indicative of diseases of the urinary system as well as other organ systems, including liver function, acid base status and carbohydrate metabolism.

### INSIGHT® URINALYSIS READER

- Semi-automatic urine test strip reader with touch screen & printer to be used with InSight Urine Test Strips
- Simply dip the strip into the urine sample and place it on the strip loading plate
- Results displayed in 60 seconds



**Use InSight Reader in conjunction with InSight MS Strips**

### INSIGHT® URINE TEST STRIPS

- Quantitative Test for the Detection of Proteinuria in Animals
- Urine Protein/Creatinine Ratio (UPC)
- Screens patients for the early signs of renal disease and monitor thereafter

MS-11 Urine Strips Test for:

- Urobilinogen
- Blood
- Glucose
- pH
- Bilirubin
- Protein
- Ketone
- Creatinine
- Specific Gravity
- Nitrate
- Leucocytes

with urine protein/creatinine ratio



MS-2 Urine Strips Test for:

- Microalbumin
- Creatinine

with urine microalbumin/creatinine ratio



**UPC ratio results in less than 2 minutes**



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT  
LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**