



qLabs[®] Veterinary Coagulation Analyser

Calibrated and validated specifically
for use in dogs and cats



COAGULATION SCREEN.
ONE DROP WHOLE BLOOD.

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.



qLabs[®] Veterinary Coagulation Analyser

In-Practice Point of Care Coagulation Testing

Calibrated and validated specifically for use in dogs and cats

KEY FEATURES

- Requires 1 drop of non-anticoagulated whole blood (10µl)
- PT/aPTT tested simultaneously on 1 test strip
- Room temperature test strip storage
- Fill test strip directly from syringe
- QC material available
- No need to transfer sample into citrate tube
- 'Test count up' for immediate assessment

SIMPLE TEST PROCEDURE

1. Insert test strip
2. Apply a drop of blood
3. Results displayed on LCD screen

HELPING DIAGNOSE CONDITIONS IN RELATION TO:

- Pre-operative checks for at risk animals
- Toxins and systemic diseases
- Hereditary bleeding disorders
- Excessive bleeding
- Disseminated Intravascular Coagulopathy (DIC)

CARTRIDGE TYPES

- PT/aPTT combo cartridge (12)

KIT CONTENTS

- qLabs[®] Coagulation Analyser
- Power Adapter
- Lancets (10)
- Lancing Device
- User Manual
- Carry Case



2 Year Warranty

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**