

www.woodleyequipment.com



:: RAPID DIAGNOSTIC TEST KITS ::



The Woodley InSight™ range of veterinary rapid diagnostic tests are simple and easy to use and provide accurate results in just 10 minutes. Based on sandwich lateral flow immunochromatographic assays, the test kits contain 10 devices per box and are conveniently stored at room temperature (2^oC - 30^oC)

InSight™ Feline Calicivirus (FCV)

InSight™ FIV- FeLV Combi

InSight™ FeLV

InSight™ FIV

InSight™ Canine Parvovirus

InSight™ Giardia

InSight™ Leishmania Ab

InSight™ Heartworm

Tel: +44 (0)1204 669033 (Option 1)

Fax: +44(0)1204 669034

Email: sales@woodleyequipment.com

Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK

InSight™

RAPID TESTS

:: THE WOODLEY INSIGHT RANGE ::

InSight™ Feline Calicivirus (FCV)

Rapid diagnostic test kit for the detection of Feline Calicivirus (FCV) in saliva, nasal and eye discharge of sick cats

InSight™ FIV - FeLV Combi

Rapid diagnostic test kit for the detection of feline immunodeficiency antibody and feline leukaemia antigen in cat serum, plasma or whole blood

InSight™ FIV

Rapid diagnostic test kit for the detection of feline immunodeficiency antibody in cat serum, plasma or whole blood

InSight™ FeLV

Rapid diagnostic test kit for the detection of feline leukaemia antigen in cat serum, plasma or whole blood

InSight™ Canine Parvovirus

Rapid diagnostic test kit for the detection of canine Parvovirus in dog faeces or vomit

InSight™ Giardia

Rapid diagnostic test kit for the detection of Giardia antigen in animal faeces


InSight™ Leishmania Ab

Rapid diagnostic test kit for the detection of Leishmania antibody (LSG Ab) in dog's serum or plasma


InSight™ Heartworm

Rapid diagnostic test kit for the detection of heartworm antigen in dog serum, plasma or whole blood

Simple Test Procedure



Interpretation of Results



Plus.....more tests in development

Tel: +44 (0)1204 669033 (Option 1)

Fax: +44(0)1204 669034

Email: sales@woodleyequipment.com

Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK