

InSight V3[®] PLUS

Veterinary Haematology Analyser



3 PART DIFFERENTIAL.
WITH EOSINOPHIL FLAG.

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

FEATURES

- Quick and easy to use
- Small sample volume - 8.5µl whole blood
- 60 second analysis time
- 3-part differential, 21 parameters with 3 histograms and eosinophil flag
- Species specific reference ranges
- Impedance technology
- Stores up to 100,000 results
- Fully automated
- Touch screen
- Low maintenance
- Results transmitted via InSight p-LINK to Practice Management Systems



INTEGRATES WITH PRACTICE MANAGEMENT SYSTEMS

SPECIFICATIONS

Principles

- Uses Coulter electrical impedance and colourimetry methods to test the WBC, RBC, PLT and HGB parameters
- 3 WBC differentials, with histograms and eosinophil flag

Parameters:

WBC	GRAN#	MCHC	PCT
LYM%	RBC	RDW_CV	P_LCR
MID%	HGB	RDW_SD	P_LCC
GRAN%	HCT	PLT	
LYM#	MCV	MPV	
MID#	MCH	PDW	

Species:

Dog, Cat, Horse, Pig, Cow, Buffalo, Rabbit, Monkey, Rat, Mouse, Sheep, Goat, Camel

User Interface

- Easy to use, menu driven user interface
- Screen Type - 10.4" capacitive colour touch screen LCD

Connection Interfaces

- LAN interface (1 x RS-232 Port)
- Printer Options - Built-in thermal printer or connect to external inkjet and laser printers
- 3 x USB interfaces

Graphic Parameters

- Displays results of 21 parameters
- Displays WBC, RBC and PLT histograms



Software Languages

English, Chinese, Spanish, French, Russian, Kazakh, Portuguese, Turkish, Italian, Czech

Dimensions

- Width 30.5cm (12")
- Depth 38cm (15")
- Height 39.5cm (15.5")
- Weight - 18kg (39.7lbs)
- Power Supply - 100-240V AC

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**