

# **InSight V-IA<sup>®</sup> PLUS**

## **Veterinary Immunoassay Analyser**

**An easy to use Immunoassay Analyser providing accurate and reliable results in 3-15 minutes.**



**LOW COST,  
REFERENCE LAB QUALITY,  
IMMUNOASSAY TESTING.**



**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**

# InSight V-IA<sup>PLUS</sup> Veterinary Immunoassay Analyser

An easy to use Immunoassay Analyser providing accurate and reliable results in 3-15 minutes

## FEATURES

- Immunofluorescence technology
- Testing time 3-15 minutes (depending on analyte)
- Easy to use
- Accurate and reliable
- Multi-language operation system
- Wi-Fi connectivity
- Reagents stored at room temperature
- Touch screen
- Maintenance free
- Results transmitted via InSight p-LINK to Practice Management Systems



### Software Languages

English, Chinese, Spanish, French, Russian, Portuguese, Turkish, Italian, Polish, German

## SPECIFICATIONS

**Dimensions:** 213 (W) x 243 (D) x 195 (H) mm  
**Weight:** 4kg  
**Power Adapter:** 100~ 240V, 50~ 60Hz

### Tests In Development:

- fHbA1C
- HP (Helicobacter Pylori)
- CPSE
- Fructosamine
- Equine ACTH
- Canine Relaxin
- Giardia Ag
- Rotavirus Ag

### Canine Tests Available:

- cPL
- T4
- TSH
- Cortisol
- Progesterone
- NT-proBNP/Troponin I (combined test)
- Cystatin C (SDMA alternative)
- Alpha-fetoprotein
- cCRP
- D-Dimer
- HbA1C
- NT-proBNP
- Troponin I
- **Core Antibody:** Distemper, Parvovirus, Distemper/Parvovirus/Infectious Hepatitis Virus combined test, Rabies Virus, Toxoplasma Gondii
- **Infectious Disease:** Coronavirus, Distemper, Parvovirus, Parvovirus/Coronavirus combined test

### Feline Tests Available:

- fPL
- fTT4
- TSH
- Cortisol
- NT-proBNP/Troponin I (combined test)
- Cystatin C (SDMA alternative)
- **Core Antibody:** Calicivirus, Herpes, Panleukopenia Virus, Herpes/Calicivirus/Panleukopenia Virus combined test, Rabies Virus, Toxoplasma Gondii
- **Infectious Disease:** Coronavirus, Calicivirus, Herpes, Panleukopenia Virus,
- D-Dimer
- fSAA
- NT-proBNP
- Troponin I

### Equine Tests Available:

- eSAA
- Progesterone
- Foal IgG



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT  
LABORATORY DIAGNOSTICS"**

**WOODLEY  
EQUIPMENT COMPANY LTD.**