

# RAPIDPoint® 500e

## Blood Gas System



**CRITICAL CARE.**  
**TRUSTED TECHNOLOGY.**



Advanced  
Partner

Laboratory Diagnostics  
Point of Care Solutions

**SIEMENS**  
**Healthineers**

**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY**  
EQUIPMENT COMPANY LTD.

# RAPIDPoint® 500e

## Blood Gas System

Siemens Healthineers RAPIDPoint 500e Blood Gas System is designed to meet the challenges of Point of Care testing by providing fast, accurate and comprehensive test results.

### FEATURES

- Accurate, lab-quality test results in 60 seconds
- Easy to use, reliable, maintenance-free operation with cartridge based system, automatic calibration and QC
- Robust accuracy for sample to sample by leveraging the Siemens Healthineers designed Integri-sense™ technology
- Onboard software defences provide reassuring level of protection against cyber attacks and confidential data theft
- Simple touch screen - simply press start after inserting a sample
- Preset custom panels
- Hands-free automatic sampling with universal sample port
- Integrated barcode scanner



**Ideal for sites with the larger workload**

### RAPIDPoint 500e Tests

**Blood Gas:** pH, pCO<sub>2</sub>, pO<sub>2</sub>  
**Electrolytes:** Na<sup>+</sup>, K<sup>+</sup>, Ca<sup>++</sup>, Cl<sup>-</sup>  
**Metabolites:** Glu, Lac  
**Co-oximetry:** nBili, tHb, sO<sub>2</sub>, O<sub>2</sub>Hb, HHb, COHb, MetHb  
**Sample Volume:** 100µl

### Weight & Dimensions

Weight: 16.55kg  
Width: 30cm  
Depth: 42cm  
Height: 55cm  
Touchscreen: 21.1x15.8 cm

### Communications

- LIS
- Wireless
- Dual-port transmission via Ethernet and serial port

### Data Security

Depend on the RAPIDPoint 500e System to raise the bar in blood gas IT security with the latest defences to guard confidential patient data coupled with the leading built-in technology to protect from external cybersecurity threats. The Windows® 10 operating system has been equipped with McAfee embedded anti-malware, which blocks unknown and unauthorised programs from running. Other key security advantages include a two-step authentication process, an encrypted password requirement for patient data transfer and the inclusion of a firewall to block hacking attempts.

### Integri-sense™ Technology

Integri-sense™ Technology is a comprehensive series of analyser functional checks and flagging mechanisms designed to deliver accurate test results. It performs frequent quality and blood integrity checks before, during and after every patient sample. Three levels of independent automatic quality control, multiple calibration routines and advanced software algorithms combine to enable the RAPIDPoint 500e Analyser to be ready to generate reliable and clinically actionable test results each and every time.

### EASY OPERATION



1. Scan the sample using the integrated barcode reader



2. Insert the sample



3. Analyse and receive results in 60 seconds

Advanced Partner

Laboratory Diagnostics  
Point of Care Solutions

**SIEMENS Healthineers**

+44 (0) 1204 669033 (Option 1)

sales@woodleyequipment.com

www.woodleyequipment.com

Old Station Park Buildings  
St. John Street  
Horwich  
Bolton  
Lancashire  
BL6 7NY  
UK

**"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"**

**WOODLEY EQUIPMENT COMPANY LTD.**