

THE SPOTCHEM EL

:: ELECTROLYTE ANALYSER WITH ION-SENSITIVE ELECTRODES ::



- Easy calibration
- Small sample volume
- Convenient twin-pipette
- Simultaneously tests electrolytes Na, K and Cl
- Highly reliable electrodes
- Low maintenance
- Test results in 60 seconds

**FAST, ACCURATE AND SIMPLE
ELECTROLYTE ANALYSIS OF BODY FLUIDS
ESSENTIAL IN AN EMERGENCY**

3 Simple Steps



Set the exclusive E-Plate



Draw the sample and reference solution with the Twin-pipette. Set pipette and apply sample to E-plate



Remove pipette. Measurement starts automatically. Results in one minute

Sample	Whole blood, plasma, serum and urine
Electrolytes	Na, K and Cl
Measurement Principle	Potentiometry, by ion-selective electrodes
Measurement Time	60 seconds
Sample Volume	22µl
Data Storage Volume	50 results
Printer	58mm width built-in thermal printer
External Output	RS-232C (serial)
Power Supply	AC100-240V 50/60Hz
Power Consumption	40VA
Weight	1.5kg
Dimensions	138h x 135w x 225d mm



Tel: +44(0)1204 669033
Fax: +44(0)1204 669034
Email: sales@woodleyequipment.com
Web: www.woodleyequipment.com

WOODLEY
EQUIPMENT COMPANY LTD.

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HORWICH, BOLTON, LANCASHIRE, BL6 7NY, UK